

BCI standards

Christoph Guger
g.tec, Austria, Spain, USA





MOZART



Musical
Empress
Elisabeth



Emperor's castle



Research Projects

H2020 SME project: **recoveriX** - motor recovery after stroke



H2020 SME project: **ComaWare** – coma assessment and communication



H2020 Eurostars project: **ComAlert** – coma prediction

H2020 Eurostars project: **RapidsMaps** – high gamma mapping

EC project: **Neurographene** – development of Graphene electrodes



EC project: **ReNaChip** - Rehabilitation of a discrete sensory motor learning function



EC project: **Sm4all** – Smart Home for all



EC project: **RGS** – Rehabilitation Gaming System faster recovery from stroke



EC project: **BrainAble** - BCI with VR and social networks



EC project: **Decoder** - BCI for locked in patients



EC project: **CSI** - Central Nervous System Imaging

EC project: **BETTER** BCI for Stroke rehabilitation and rehabilitation robots



EC project: **VERE** – Virtual Embodiment Real Embodiment



EC project: **ALIAS** – Adaptable Ambient Living Assistant



EC project: **BACKHOME** – BCIs for end users




EC project: **DENECOR**



EC project: **High Profile**



Current g.tec BCI applications



Changes of brain electrical activity and mental strategies

- Slow cortical potentials (anticipation tasks)
DC-derivation, artifact problem, difficult strategy, feedback method
- Steady-State Evoked potentials (SSVEP, SSSEP)
Flickering light with specific frequency
- Event-related, non-phase-locked changes of oscillatory activity
ERD/ERS (motor imagery tasks)
Changes of mu-rhythm, alpha activity and beta activity over sensorimotor areas;
imagination of hand-, foot-, tongue- movements
- Evoked potentials (focus on attention task)
Thalamic gating, various methods of stimulation (visual, tactile, electrical, auditory, ...),
P300
- Code based evoked potentials (focus on a code)

Problem 1: Different sensors

Non-invasive versus invasive (different regulations, FDA, CE)

Depth electrodes, grids versus EEG electrodes

Active or passive electrodes

Gel or dry EEG electrodes

Different number of channels (P300: 8, SSVEP: 8, motor imagery: 64, invasive mapping: 256)

Different sampling frequency (Spikes: 40 kHz, ECoG: 1-4 kHz, EEG: 256 Hz)

Different platforms (Windows, Linux, Android,...)



- A** g.MOBilab+ ☐ EEG version
☐ multi-purpose version

☐ g.GAMMACap



☐ Consumables

- Electrodes ☐ active
☐ passive
☐ dry



- B** g.USBamp ☐ 16 channels
☐ 32 channels
☐ 48 channels
☐ 64 channels



- C** g.HIamp ☐ 80 channels
☐ 144 channels
☐ 256 channels



- D** g.Nautilus ☐ 8 channels
☐ 16 channels
☐ 32 channels
☐ active-gel
☐ active-dry



- ☐ Eyetracker (screen-mounted)
☐ Eyetracker (head-mounted)



Sensors

- ☐ GSR
☐ Blood pressure
☐ Respiration
☐ Temperature
☐ Pulse
☐ Acceleration
☐ More sensors

☐ Electrical stimulator




- ☐ g.STIMbox
☐ g.SSVEPbox



Recommended setup for a fully equipped BCI lab

Room dimension: >5x6m




g.tec medical engineering GmbH, office@gtec.at, www.gtec.at, tel. +43 7251 22240, fax +43 7251 22240 39 g-tec

- ☐ Impedance tester
☐ Signal generator



☐ VR-system



Software

- | | |
|---|---|
| <input type="checkbox"/> g.Recorder | Examples |
| <input type="checkbox"/> g.BSanalyze | <input type="checkbox"/> P300 |
| <input type="checkbox"/> g.Hlsys (Simulink) | <input type="checkbox"/> Motor imagery |
| <input type="checkbox"/> g.Hlsys (LabVIEW) | <input type="checkbox"/> SSVEP |
| <input type="checkbox"/> intendiX | <input type="checkbox"/> CSP/RehaBCI |
| <input type="checkbox"/> BCI2000 | <input type="checkbox"/> Vibro-tactile P300 |
| <input type="checkbox"/> C/MATLAB/Linux API | <input type="checkbox"/> Ping-Pong game |
| <input type="checkbox"/> g.RTanalyze | <input type="checkbox"/> Hyperscanning |
| <input type="checkbox"/> g.PHYSIOobserver | <input type="checkbox"/> Hybrid BCI |
| <input type="checkbox"/> g.UOPInterface | <input type="checkbox"/> EMG/EOG control |
| <input type="checkbox"/> mindBEAGLE | Training material |
| <input type="checkbox"/> cortiQ | <input type="checkbox"/> EEG lecture |
| <input type="checkbox"/> Eyetracker interface | <input type="checkbox"/> EP lecture |
| <input type="checkbox"/> g.NEEDaccess | <input type="checkbox"/> BCI lecture |

Request an offer for a complete lab
☐ minimum ☐ basic ☐ advanced

Check products of interest and request an offer per fax (+43 7251 22240 39) or e-mail (office@gtec.at) !

How do interface with different amplifiers

g.NEEDaccess service to interface all amplifiers with one common interface

A. C++ Application Program Interface (API)

- integrate amplifiers into own software under Windows and Linux

B. MATLAB API


- integrate amplifiers into MATLAB data acquisition and analysis programs
- access all toolboxes (Signal Processing, Neural Networks,...)
- access user written M-files

C. Simulink Highspeed on-line Pro

- amplifier device driver block
- copy the block into Simulink (S-functions) and paradigm
- just exchange the amplifier processing blocks

D. LabView

- amplifier device driver block
- use standard LabView blocks



All three options give full access to hardware

- bandpass, notch settings
- sampling frequency
- impedance check
- synchronization with digital inputs
- direct integration of other devices

Problem 2: Event Timing

Real-time data stream synced with external devices

visual P300: 1 ms resolution

vibro-tactile P300: 1 ms resolution


....

Control of external devices

UDP


digital outputs

intendiX – face speller





**Perfect timing between flash and EEG
needed (1ms resolution) -> signal
processing must be perfectly linked to
visual presentation.**

- 17 subjects
- 4 tasks
 - Standard flashing
 - „Einstein“ (black/white)
 - „Einstein“ (color)
 - Face speller



Everybody achieved 100 % accuracy

BCI interface with video overlay





Code stimulation must be perfectly timed -> refresh rate of the monitor is essential

Video-Server / Video-Client / BCI-Overlay / LEDs
(visual stimulation and video feedback)

11 subjects participated

LED stimulation: 91,36 %
Screen stimulation: 91,36 %
Code based screen: 98.18 %.




Average time to complete the tasks
222.57 s (code based BCI),
437.43 s (frequency LED)
573.43 s (frequency screen).



Submitted to Frontiers, 2013

Robot control in VERE [video](http://www.gtec.at)
www.gtec.at

Embodiment Station



Gerwin Schalk, Wadsworth Center

Stimulate the body and observe effects in the brain -

-> real-time functional mapping

Stimulate the brain and observe effects on the body

-> real-time sensing

Aim: Build a functional cortical atlas.

[cortiQ - Clinical software for electrocorticographic real-time functional mapping of the eloquent cortex.](#)

Prueckl R, Kapeller C, Potes C, Korostenskaja M, Schalk G, Lee KH, Guger C.

Conf Proc IEEE Eng Med Biol Soc. 2013 Jul;2013:6365-8. doi: 10.1109/EMBC.2013.6611010.


Problem 3: Many controllable elements

Smart home needs many controls for domotic devices

BCI has to understand which controls are necessary

BCI must be updated to changes of the environment

Brain Painting Application




Designed by Adi Hösle

www.gtec.at


ACTOR protocol

- The BCI speaks with ACTOR protocol with the avatars and robotic systems
- **XML files** are loaded **at startup** (from disk or from a text string, which is received over the network)
- **Updates** of the XML files can be received at runtime over UDP: **Modify the contents of the BCI at runtime**, e.g. to achieve context awareness



The XML based BCI Standard

- Each **interaction** (e.g. the letter 'A') is related to an instruction string
- If 'A' is selected: the **string is sent over UDP**:
A **remote device or applications can be controlled with the BCI** just by knowing an IP address and the instruction string



Send instruction 'A' to a remote device or application (e.g. Twitter)

The Virtual Reality Apartment




The XVR 3D-software allows to create highly immersive virtual scenarios and to control them via a simple and flexible interface. Physiological parameters of the subject (such as EEG, EMG, ECG, EOG, GSR and respiration) can be recorded and processed in real-time.



award-winning technology.

contact: g.tec medical engineering GmbH, Sierningstr. 14, 4521 Schiedlberg, AUSTRIA, phone: +43 72

In a recent research project a BCI (Brain-Computer Interface) is used to navigate in a virtual Smart Home and to control various environmental parameters (e.g. music, TV, temperature, light, doors or games). The required biosignals (EEG) are transmitted wirelessly to a BCI computer where data processing and analysis takes place. The output commands generated by the BCI controls the VR via a UDP interface.



Problem 4: Avatar/robot control

The person is seeing the environment through the avatar or robotic system

We need BCI controls for controlling the avatar/robot

The BCI system has to send control command to external system

Interfacing with rehabilitation devices

Screen overlay control interface - SOCI

Video: [Overlay BCI](#)

World of Warcraft



4 controls:

Turn left, right, move forward, perform action like grasping objects, attacking other objects

60 Hz LCD display with 15, 12, 10 and 8.75 Hz.

BCI overlay based on OpenGL –
can be used with any graphics application



Humanoid Robots for Physical Embodiment



[BCI control to grasp a Coke](#)



[EuroNews](#)



[DigInfo](#)

Abderrahmane Kheddar, CNRS

Problem 5: Highly immersive feedback



Stroke rehabilitation with BCI system

Motor imagery BCI controls FES and avatar

BCI accuracy is objective parameter for training

Motor recovery controlled with 9-hole PEG test



Problem 6: Performance standards






BCI performance comparison

**Important that some systems are not called
BCI systems**


	Motor imagery	cVEP	P300 speller [Guger 2016]
Grand average accuracy	87 %	98 %	100 %
Training time	30 min	5 min	5 min
Number of electrodes	32	8	8
Random classification accuracy	$\frac{1}{2}$	$\frac{1}{4}$	$\frac{1}{36}$
Decision time for selection	6 sec	3 sec	About 45 sec with 15 flashes
Location	Motor cortex	Visual cortex	Central line and visual cortex

Hand Movement Task

- ➞ Movement: Online classification accuracy ([video](#))
- ➞ Imagination: Online classification accuracy ([video](#))



BCI-ECoG setup




Classification accuracy

Problem 7: Calibration



SSVEP group study accuracy

Accuracy (%)	Number of subjects performing at specified accuracy				Percentage of people after training
	Run 1	Run 2	Run 3	Run 4	
100	22	25	27	27	50.9
90-99	14	19	BCIs need calibration phase to achieve high performance		
80-89	7	4			
70-79	2	1			
60-69	1	2	1	1	1.9
50-59	4	1	0	0	0.0
40-49	3	0	1	0	0.0
0-39	0	1	0	0	0.0
Mean Accuracy	87.9	92.9	95.0	95.5	
	N=53	N=53 with 14 new	N=53 with 7 new	N=53 with 2 new	



Poor performance in SSVEP BCIs: Are worse subjects just slower?



frontiers
IN NEUROPROSTHETICS

How many people could use an SSVEP BCI?, Christoph Guger, Brendan Z Allison, Bernhard Grosswindhager, Robert Prückl, Christoph Hintermüller, Christoph, Kapeller, Markus Bruckner, Gunther Krausz and Guenter Edlinger, Frontiers in Neuroprosthetics, 2012.



Problem 8: Good looking BCI users

